



THE 2HR AQUARIST



APT
handling data

plant supplements

APT1 APT ZERO

	Q1 2022				
Size	100ml	200ml	300ml	500ml	1000ml
Code	APT1-100	APT1-200	APT1-300	APT1-500	APT1-1000
EAN	NA	8881300100049	8881300100148	8881300100056	8881300100063
N Weight	NA	225g	356g	562g	1129g
G Weight	NA	255g	398g	604g	1195g
Dimension	NA	50x50x170mm	59x59x200mm	61x61x210mm	82x82x222mm

APT3 APT COMPLETE

	Q1 2022				
Size	100ml	200ml	300ml	500ml	1000ml
Code	APT3-100	APT3-200	APT3-300	APT3-500	APT3-1000
EAN	NA	8881300100018	8881300100155	8881300100025	8881300100032
N Weight	NA	225g	358g	562g	1129g
G Weight	NA	255g	400g	604g	1195g
Dimension	NA	50x50x170mm	59x59x200mm	61x61x210mm	82x82x222mm

APT_e ESTIMATIVE INDEX

	Q1 2022				
Size	100ml	200ml	300ml	500ml	1000ml
Code	APTe-100	APTe-200	APTe-300	APTe-500	APTe-1000
EAN	NA	8881300100070	8881300100162	8881300100087	8881300100094
N Weight	NA	225g	378g	562g	1129g
G Weight	NA	255g	420g	604g	1195g
Dimension	NA	50x50x170mm	59x59x200mm	61x61x210mm	82x82x222mm

tank conditioners

APT_p APT PURE

	Q1 2022				
Size	100ml	200ml*	300ml	500ml	NA
Code	APT _p -100	APT _p -200	APT _p -300	APT _p -500	NA
EAN	NA	8881300100117	8881300100650	8881300100100	NA
N Weight	NA	304g	370g	598g	NA
G Weight	NA	334g	412g	640g	NA
Dimension	NA	50x50x170mm	59x59x200mm	61x61x210mm	NA

APT_f APT FIX

	Q1 2022				
Size	100ml	200ml#	300ml^	500ml#	NA
Code	APT _f -100	APT _f -200	APT _f -300	APT _f -500	NA
EAN	NA	8881300100667	8881300100131	8881300100223	NA
N Weight	NA	201g	320g	499g	NA
G Weight	NA	232g	363g	542g	NA
Dimension	NA	50x50x170mm	59x59x200mm	61x61x210mm	NA

Notes

- ^APT Fix 300ml Standard comes with dosing pipette in box, some variants with screw-cap instead of pump.
- #APT Fix 200ml and 500ml come with attached pipette.
- *Regular 200ml bottle filled to 250ml, with sticker "25% extra- 250ml"

Carton Box Size: 41cm x 21cm x 26cm
 Each carton box weighs 500g and fit:
 10 x 1000ml
 20 x 300ml
 20 x 500ml
 40 x 200ml
 40 x 100ml